



## Thermostable Reverse Transcriptase [M-MLV, RNase H-]

### User's Instruction

**Product Concentration:** 200 U/μl

**Storage:** -20°C, 2 years

#### Description

Thermostable Reverse Transcriptase is a recombinant M-MLV reverse transcriptase with reduced RNase H activity and increased thermostability. The enzyme is active up to 55°C. It provides higher specificity, higher yield and more full-length cDNA products.

- Increased thermostability for more full-length cDNA products.
- Deficient RNase H activity to reduce RNA template degradation during the first-strand cDNA synthesis.
- cDNA up to 15 kb.

#### Unit Definition

One activity unit (U) refers to the amount of M-MLV reverse transcriptase when catalyzes the incorporation of 1 nmol of dTTP into materials in 10 min at 37°C using oligo (dT) primed poly (A) as a template.

#### 5 X RT Buffer

250 mM KCl; 15 mM MgCl<sub>2</sub>; 10 mM DTT; 100 mM Tris-HCl pH 8.4.

#### First Strand cDNA synthesis (20 μl reaction volume)

1. Add components according to the below table:

Components	Volume
Total RNA/mRNA	50 ng-5 μg/5-500 ng
Oligo(dT) <sub>18</sub> (0.5 μg /μl)	1 μl
<b>Or</b> random Primer (0.1 μg/μl)	1 μl
<b>Or</b> GSP (Gene Specific Primer)	2 pmol
dNTP Mix, 10 mM each	1 μl
5 X RT Buffer	4 μl
Ribonuclease Inhibitor (40 units/μl)	0.5 μl
Thermostable H <sup>-</sup> RTase	1 μl
RNase free H <sub>2</sub> O to final volume	20 μl

**Optional** (if RNA template is GC-rich or is known to contain secondary structures). Suggest to mix RNA /Primer/RNase free H<sub>2</sub>O gently and briefly centrifuge, incubate at 65°C for 5 min, chill on ice and briefly centrifuge, then place the tube on ice. Add other components and continue.

2. Mix well gently



If Oligo(dT)<sub>18</sub> or gene specific primer (GSP) are used, incubate at 50°C for 30-50 min.

If Random Primer is used, incubate 10 min at 25°C followed by 30-50 min at 50°C.

3. Terminate the reaction by heating at 70°C for 15 min.

The reverse transcription reaction product can be directly used in PCR or stored at -20°C.

## RT-PCR

Use 2-4 µl of the reaction mix to perform PCR in 50 µl volume.

**PCR mixture set up** (for 50 µl reaction volume)

Components	Volume	Final Concentration
cDNA Template	2-4 µl	as required
Forward Primer (10 µM)	1 µl	0.2 µM each
Reverse Primer (10 µM)	1 µl	0.2 µM each
10X Taq Buffer (contains Mg <sup>2+</sup> )	5 µl	1×
2.5 mM dNTPs	4 µl	0.2 mM
Taq DNA Polymerase	0.5 µl	2.5 units
ddH <sub>2</sub> O to final volume	50 µl	Not applicable

## PCR Condition

94°C	2-5 min	} 30-40 cycles
94°C	30 sec	
50-60°C	30 sec	
72°C	1-2 kb/min	
72°C	5-10 min	