

NEXTflex™ Msp1 Restriction Enzyme

产品特点:

- Optimized for bisulfite sequencing

产品描述:

NEXTflex™ Msp1 Restriction Enzyme is a restriction enzyme that leads to selective amplification of CpG regions, resulting in focused sequence depth.

The NEXTflex™ Msp1 Restriction Enzyme is designed to be used with the NEXTflex Bisulfite-Seq Kit for reduced representation. Unlike other commercially available MSP1 restriction enzymes, the NEXTflex™ Msp 1 Restriction Enzyme has been optimized for bisulfite sequencing.

Catalog#	Product Name	Quantity
511921	NEXTflex™ Msp1 Restriction Enzyme	8 rxns
511922	NEXTflex™ Msp1 Restriction Enzyme	48 rxns