

## AIR™ DNA Fragmentation Kit

### 产品特点

- Optimized to physically shear uniform fragments of genomic DNA
- Allows smaller sample sizes and reduces shearing time
- Offers 30% > DNA recovery than competitor's fragmentation protocol

### 产品描述

The AIR™ DNA Fragmentation Kit is designed to shear genomic DNA into 150 to 800 bp fragments which are ideal for second generation sequencing technologies. Nebulization repeatedly forces input DNA through a small hole in the unit producing reproducibly random fragments. These fragments are double stranded and contain 3' overhangs, 5' overhangs and blunt ends in a heterogeneous mix.

Catalog#	Product Name	Quantity	Price
5135-01	AIR™ DNA Fragmentation Kit	10 rxns & nebulizers	Inquire
5135-02	AIR™ DNA Fragmentation Kit	40 rxns & nebulizers	Inquire